

SWISS GROUP TEAM



AMBULANCE CONVERSION

SPECIFICATIONS AND DETAILS

The following medical equipment and items will be fitted into a Level C conversion ambulance

MEDICAL EQUIPMENT:

Main wheeled transporter

Main stretcher with safety belt.

The stretcher is mounted on the floor with rail and locking system. Tubular Aluminium construction provides durability & strength. This cot adjusts to six height levels. Adjustable backrest angle from 0° to 75° allows patient to be positioned during transportation. Two patient restraints ensure maximum safety of the patient during transport. One of the lightest lift-in cot in the industry, easily installed in ambulance.



| | |
|----------------------------|----------|
| Maximum Height | : 85 cm |
| Minimum Height | : 32 cm |
| Maximum Length | : 193 cm |
| Overall width | : 58 cm |
| Weight (without mattress): | 29 kg |
| Load limit | : 180 kg |

Folding Stretcher

Aluminium Alloy Folding Stretcher

Stretcher Opened: Approx. 2000(L) x 550(W) x 150(H) mm

Tubes: Approx. Dia. 30 mm x 1.5 mm

Carrying Capacity: Approx. 150 kg



Oxygen Cylinder (plus bracket) with regulator & mask



- Flow control range: 1-10l / min
- Input pressure: <15mpa
- Output pressure: 0. 2~0. 3mpa
- The control pressure of safety valve: 0. 35 0. 05mpa
- 48 cu ft Oxygen Cylinder

Single head stethoscope

Brass chrome-plated binaurals
56cm y-tube
Anti chill ring
Membrane 44 mm diameter
Brass chrome-plated chest piece



INTERIOR:

Storage area for medication

This is a typical interior as designed by us with ample storage space for medication. This can of course be adapted to suit your individual needs and requirements.



Partition between driver compartment and patient area

Partition window comprising two large sliding windows with safety glass and roller blind in driver's compartment

Patient compartment is completely insulated (ceiling, side, floor and back doors)


Patient compartment electrical system separated from track

Patient compartment has sliding windows for easy access and is 2/3rds shaded

Electric socket points (2 internally and 1 externally)

IV hooks

Cushioned bench seating for up to 4 persons with storage space



Roof ventilator

Well designed ceiling mounted one way extractor which eliminates air and odours from the patient compartment and disposes them outside the ambulance.



Examination lamp & portable search light

The examination lamp is operated by switch and is mounted flush to ceiling. The search light is 12v and is portable.



One Paramedics chair at the head of the patient compartment with 3 point seat belt

High quality flooring material

First aid kit complete:



- 2 PIC 3 Box with:
- 3 towellettes w/ disinfectants
- 2 towellettes w/ ammonium hydroxide 4%
- 3 hygienic soap pouches
- 2 adhesive plaster rolls M.5 x 1.25cm
- 1 pack of gauze bandages
- 5 packs of sterile gauze cm. 10x10 12 folds
- 5 packs of 25 pcs sterile gauzes cm 10x10
- 5 cotton wool packs GR.50
- 6 safety pins
- 4 triangular bandages TNT cm.96 x 96 x 136
- 1 pair of scissors 10cm
- 2 sterile tweezers
- 1 disposable sterile scalpel
- 1 haemostatic lacet
- 8 hypodermic syringes cc. 5 w/ towelettes
- 1 plastic kidney dish

4 eco splints 22cm
1 pack of 20 adhesive plasters

Air conditioning system in patient's compartment

Electronic siren with public address system of 100 watt output power

6 rotary roof beacons mounted on the roof of patient's compartment (2 at the front, 2 at the rear and one on each side)

100 watt speakers embedded on the roof at the front of the vehicle

Two fire extinguishers (one in the driver's cabin and one in patient's compartment)

Nurse chair adjacent to medical equipment

Small waste basket – sharps container

Ambulance blanket

Dustpot

Trauma Bag

Hazard Warning Triangle

Vacuum Splint Set

Airway Set

Standard ART Work:

“Ambulance” in English in red colour and Red Crescent shown on each side of the vehicle, one at the front, one at the rear and another on the roof.

Artwork may be altered to your requirements but additional charges may apply



***ADDITIONAL ITEMS MAY BE ADDED AT AN EXTRA CHARGE
PLEASE ASK FOR A QUOTATION***



www.swisstoyota.com

info@swisstoyota.com